

## **UCMSC Product Information**

	General Details		
1	Cell Name	Umbilical Cord-derived Mesenchymal Stem Cell	
	Description	Human Mesenchymal Stem Cells were isolated from umbilical cord Wharton's jelly of a healthy woman during childbirth and manufactured using a Good Manufacturing Practice (GMP) process	
	Passage number	Passage 3 / Passage 5	
	Number of Cells	4X10 <sup>6</sup> cells/vial	
	Culture Type	Adherent	
	Culture Condition	37°C, 5% CO <sub>2</sub> incubator	
	Culture Media	αMEM + 5% HPL +3 mM Glutamax	
	Storage Condition	Liquid Nitrogen	
	Shipping Condition	1 mL of Cryovial	
	Cell morphology		
	Test Specification		
	Cell Performance Testing		
2	Cell Viability	>80%	
	Chondrogenic, Osteogenic, Adipogenic differentiation potential	Positive	
	Microbial Testing		
	Sterility	Negative	
	Mycoplasma	Negative	
	Endotoxin	<0.05 EU/mL	
	Virus Testing		
	EBV, CMV, HHV-6, HHV-7, HHC-8, HLV-1, HLV-2, HAV, HCV, HBV, HIV-1, HIV-2, X-MuLV, retrovirus	Not detected	
	Cell Surface Marker Testing		
	Positive markers (CD90, CD105, CD73, CD29, CD49e, CD44, CD166)	≥95%	
	Negative markers (CD14, CD19, CD34, CD45, HLA-DR)	≤2%	



	Protocol/Guidelines	
3		<ul> <li><u>Cell Thawing and Seeding</u></li> <li>1. Thaw and seed cells at a density of 1,500 cells/cm<sup>2</sup></li> <li>2. Expand cell cultures for 5-6 days until they have reach &gt;80% confluency</li> <li>3. On Day 1, perform complete media change (100%). On Day 4 after seeding, perform a 50% media change</li> </ul>
	Cell Culture	<ul> <li><u>Cell Expansion</u></li> <li>1. Subculture the cells at a density of 1,500 cells/cm<sup>2</sup></li> <li>2. Expand cell cultures for 5-6 days until they have reach &gt;80% confluency</li> <li>3. On Day 3 and Day 6 after subculture, perform 50% media change</li> </ul>
		<ul> <li><u>Cell Harvest</u></li> <li>1. Centrifuge the cells at 330 x g for 7 min at 18°C</li> <li>2. Resuspend the cell pellet in CryoStor®CS2 cryopreservation medium at concentration of 4X10<sup>6</sup> cells/vial</li> <li>3. Store the vial in a freezing container overnight at -80°C before transferring it to liquid nitrogen</li> </ul>
	Additional Information	
4	UCMSCs were manufactured from donated umbilical cord materials in the GMP Unit	